


**C2A**



IMPORTANT, RETAIN FOR FUTURE REFERENCE: READ CAREFULLY

**ASSEMBLY & INSTRUCTION MANUAL**



## CAUTION

BEFORE using this Knee Scooter, read this manual and save for future reference.

**DO NOT** attempt to adjust or operate Knee Walker without carefully reading all instructions in this manual.

Check ALL parts for shipping damage before use. In the case of

shipping damage, Do NOT USE.

**Do NOT** install or use this equipment without first reading and understanding this instruction sheet.

If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional.

Attempting to install this equipment - Otherwise, injury or damage may occur. Before using the Knee Walker, you should be trained by a healthcare professional or dealer.

**Do NOT** use the Knee Scooter by yourself until you have been properly trained, you have practiced with someone, and you feel confident you can maneuver by yourself.

The user should ALWAYS consult with his/her physician or therapist to determine proper adjustment and usage.

## WARNINGS & SAFETY INSTRUCTIONS

Maximum Weight Capacity is 300 lbs. DO NOT exceed the maximum weight capacity (including basket contents) or serious injury could result.

Maximum weight capacity for the Knee Scooter basket is 10 lbs.

DO NOT use the Knee Scooter on stairs, curbs, obstacles, etc. Serious risk of fall or injury may occur.

Watch out for: cords, thrown rugs, toys on the floor, spilled water, or oil.

DO NOT use the Knee Walker if the brake system is not working properly.

DO NOT use the Knee Walker unless all parts are secure and moving parts are in good working order.

DO NOT lean forward, backward, or to either side while using.

DO NOT use the Knee Scooter unless seat lock lever is securely locked into place.

DO NOT lean forward while applying the brakes.

DO NOT perform any adjustments to the Knee Scooter while it is in use

DO NOT use the seat to transport people or objects.

While it is being transported in a vehicle.

DO NOT sit on the knee rest.

DO NOT hang anything from the Knee Scooter handles or frame at any time. This may cause the Knee Scooter to tip, resulting in bodily injury or Knee Scooter damage. Items should be placed in the basket.

DO NOT attempt to reach objects while using the Knee Scooter. Reaching for these objects will cause a change of the weight distribution and may tip over, resulting in injury or damage.

DO NOT walk backwards while using the Knee Scooter.

All wheels MUST be in contact with the floor at all times during use. This will ensure the Knee Scooter is properly balanced.

Always keep both hands on the handle grips when in use.

Only use accessories and spare parts authorized by ourself.

When turning the Knee Scooter, exercise caution and only turn while at a slow rate of movement.

The Knee Scooter is for individual use only. DO NOT attempt to push the Knee Scooter while someone is seated on the seat. Doing so may result in serious injury and damage to the Knee Scooter.

The Knee Scooter is not to be used as a transportation device. The Knee Scooter is a walking aid only.

A physical/occupational therapist should assist in the height adjustments of the Knee Scooter for maximum support and correct brake activation.

## CARE & MAINTENANCE

Clean with damp cloth and non-abrasive cleanser.

Clean wheels with warm water and mild cleanser. Dry with clean rag.


Ensure that all attached hardware and seat lever are always secured.

Inspect the Knee Scooter periodically to ensure the brakes, seat lever, steering, and wheels are functioning properly and that all parts are secure.

Replace any broken, damaged, or worn items immediately.


Ensure that the brakes are working correctly. Over time the brake cables may stretch with use, they can be adjusted by a qualified mechanic. Failure to do this may result in injury to the user.

## 1. PARTS INCLUDED



- ① Detachable front wheel x1
- ② Walker frame x1
- ③ Handlebar x1
- ④ Knee pad platform x1
- ⑤ Five-pointed star screw x3
- ⑥ Bag x1


## 2. INTRODUCTION



### 3.SETUP INSTRUCTIONS


#### Step 1: To Assemble Knee Walker Frame

- Unpack the knee walker from the carton.
- Rotate the lever clockwise to tighten the quick release clamp and pull the quick release clamp down to point towards the frame tubing to lock it securely.



#### Step 2: To Raise the Steering Tiller


- Align and straighten the steering column frame to the erect position. Slide the locking pin to the side, this allows the steering column to reach its full upright position. Release the pin to lock the steering column in place.
- Move locking pin into the groove on steering column. Then, depress the locking handle to secure the steering column.



### Step 3: To Adjust Height of Handlebars

Handlebars should be set at wrist level to keep your back straight.

- Remove star screw.
- Extend handlebars up to desired height and align holes in tiller;
- Tighten star screw securely.




### Step 4: To Adjust Height of Knee Platform

Stand on your good leg, hold on to a steady surface for support, and bend the leg that will be on the knee platform at a 90-degree angle.


Adjust the knee platform to meet the height of your bent leg. When you place your bent leg on the knee platform, your uninjured foot should comfortably reach flat on the floor, and your hips should be even. You should be in an upright position as if you were standing normally.

- Remove star screw.
- Insert the knee platform post.
- Extend the knee platform to desired height.
- Then, insert star screw and tighten knob securely.



### Step 5: To Attach Bag


Take out the cloth basket from the carton and hang the cloth basket iron ring into the hanging piece on the front support tube by hand; Check that it is tight and does not shake.



### Step6: To Adjust the Brakes


The brake tightness can be adjusted by the nuts at both ends of the brake line.

The brake tightness can be adjusted with the nuts at either end of the brake cable, or nuts on the rear wheel spring. Turning the screw clockwise will loosen the brake and turning the screw counterclockwise will tighten it.



### Step 7: Brake Operation

- 1). Squeeze both hand brakes lever at same time, the knee walker is parked. (PS: one handle brake only control one side parking system)
- 2). Permanent brake
  - a. Squeeze the hand brake lever and push down the parking brake lock, and—the Knee Walker will be locked.
  - b. Squeeze the hand brake lever again to release the parking brake lock.



### Ongoing Maintenance

With routine use, nuts and bolts on the steering column and brake assembly can become loose. To ensure safe operation of this scooter it is important to periodically inspect the knee walker for loose nuts and bolts and re-tighten as appropriate.

### To Operate the Knee Walker

- q. Place your hands on the handlebars.
- b. Place the injured leg comfortably resting on the center of the knee platform.
- c. Place the injured leg comfortably resting on the center of the knee platform.
- d. The Knee Walker can make turns by simply moving the handlebars with a minimum amount of effort. Turns should only be performed while at a slow rate of movement. To turn right, rotate handlebars to the right. To turn left, rotate handlebars to the left.
- e. With practice, you should be able to move forward without deviating to the opposite side.